

EVACUATED TUBES



13x100	13x75	Additive	Cap Color	Tube Type
5 ml	2/5 ml	Clot Activator With Gel	Yellow	SST
6 ml	4 ml	Clot Activator	Red	Serum
6 ml	3 ml	K2 EDTA	Purple	EDTA
6 ml	3 ml	K3 EDTA	Purple	EDTA
6 ml	4 ml	Lithium Heparin	Green	Heparin
5 ml	3 ml	Lithium Heparin with Gel	Green	PST
-	2/7 ml	3.2% Sodium Citrate	Blue	Coagulation
-	3 ml	Naf	Gray	Glucose
-	2/4 ml	3.8% Sodium Citrate	Black	ESR



Complied With ISO 6710; EN 14820

Expected draw volume throughout the entire shelf life

High efficiency separation gel and additive

Sterilized by radiation

Options:

Available as Plain or With reagents



100 tubes per shrinking packs

12 shrinking packs per case (1200)